




10x and cellranger

Jan Kubovčíak
with acknowledgements to Deepak Tanwar and SIB

GEMs

- way to label cells separately
- Gel Beads-in-emulsion



Single Cell
GEMs


Collect



Barcode




Pool, remove oil





Sequencing-ready
libraries


Chromium instrument



10,000+ GEMs in 4 minutes.



All captured **transcripts** from **single** cell:
identical + unique barcode



- reverse transcription
- breaking GEMs
- fragmentation
- primer ligation
- index PCR



Illumina sequencing

Sequencing output

	sample#	read type	
ETV6-RUNX1_1_S1_L001_R1_001.fastq.gz			
ETV6-RUNX1_1_S1_L001_R2_001.fastq.gz			
ETV6-RUNX1_1_S1_L001_I1_001.fastq.gz			

sample ID lane

For 10x all with
cellranger

- Dual indexing: second index in I2
- Indexes can also be added to fastq titles



ETV6-RUNX1_1_S1_L001_R1_001.fastq.gz

@SRR9264343.187191690 K00271:135:HK3J7BBXX:3:2110:12591:47964 length=26

GGCAATTCATTGGGCCCCGAAGCTC

+

AAF-FJFJJJJJJJJJJJJJJJJJJ

@SRR9264343.5519124 K00271:135:HK3J7BBXX:3:1112:2950:8295 length=26

GGCAATTCATTGGGTGAGGTTGTAGG

+

AAFFFJJJJFFJFJJJJJJJJJJ

@SRR9264343.149664941 K00271:135:HK3J7BBXX:3:2217:29579:19267 length=26

TCGGATGTCTGAAAGAGTGCTAGAAT

...

ETV6-RUNX1_1_S1_L001_R2_001.fastq.gz

@SRR9264343.187191690 K00271:135:HK3J7BBXX:3:2110:12591:47964 length=98

GTAGTAAAAGACTGGTTAATGATAACAATGCATCGTAAACCTTCAGAAGGAAGGAGAATGTTTTGTGGACCACTTTCGTTTTCTTT
TTTGC GTGTG

+

-777<<A-<FF<<A7-<FAFAFFJJJJJFFFFFJJFFJJJJFFJJJJ<FJJAFJJFJJ<A-<AA<FFF)<A--7-7--AFFJ--AFAFA7<AA-

@SRR9264343.5519124 K00271:135:HK3J7BBXX:3:1112:2950:8295 length=98

AGACTTTCTACCTGGTCATATACTCTGCAGCTGTTAGAATGTGCAAGCACTTGGGGACAGCATGAGCTTGCTGTTGTACACAGGGTA
TTTCTAGAAGC

+

--77-77--7-77<-7AFAJJJJJ7A7AFF7-A-7FFJJFA-AF77FF<-A<FFJJF-)FF7F<7FA7FAAAFJ-A--A<AAFJJ<F----7--<F77

@SRR9264343.149664941 K00271:135:HK3J7BBXX:3:2217:29579:19267 length=98

GGTCCGCAAGCAGGTGGTGCACATCCCGTCCTTCATTGTCCGCTGGCTTCCCAGAAGCGCATCGACTTCTCTCTGCGCTCTCCCTA
CGGGGGTGGCC

...

ETV6-RUNX1_1_S1_L001_I1_001.fastq.gz

@SRR9264343.187191690 K00271:135:HK3J7BBXX:3:2110:12591:47964 length=8

CCTAGACC

+

AAAFFJJJ

@SRR9264343.5519124 K00271:135:HK3J7BBXX:3:1112:2950:8295 length=8

CCTAGACC

+

AAAFFJJF

@SRR9264343.149664941 K00271:135:HK3J7BBXX:3:2217:29579:19267 length=8

CCTAGACC

...

After sequencing (pre-processing)

1. Read mapping
2. Quantification: # UMI/gene
3. Cell calling

For 10x all with
cellranger

Alternatives:

[STARSolo](#)
[Alevin](#)

cellranger references


- Human & mouse: download pre-built from 10x website
- Other organisms: custom reference with **cellranger mkref**
 - Use **ensembl.org**
- Exogenous marker genes (e.g. GFP): add sequence to both fasta and gtf

Why count UMI (and not read alignments?)

- UMI: Unique Molecular Identifier:
 - Identifies each molecule (i.e. sequence) uniquely
- Molecules from a common PCR template
-> carry the same UMI
- By counting UMI: correct for PCR replicates

Why UMIs

- Reads with **different UMIs** were derived from **different molecules** and are **biological** duplicates
- Reads with the **same UMI** originated from the **same molecule** and are **technical** duplicates – the UMIs should be collapsed to be counted as 1



Cellranger count report

ETV6-RUNX1_1

Alerts

The analysis detected ⚠️ 1 warning.

Alert	Value	Detail
⚠️ Fraction of RNA read bases with Q-score >= 30 is low	59.4%	Fraction of RNA read bases with Q-score >= 30 should be above 65%. A lower fraction might indicate poor sequencing quality. This is Read 1 for the Single Cell 3' v1 chemistry and Single Cell 5' paired end, Read 2 for the Single Cell 3' v2/v3 chemistry and Single Cell 5' R2-only)


Summary	Analysis
---------	----------

3,091
Estimated Number of Cells

68,259	1,717
Mean Reads per Cell	Median Genes per Cell




Sequencing ?	
Number of Reads	210,987,037
Number of Short Reads Skipped	0
Valid Barcodes	98.2%
Valid UMIs	100.0%
Sequencing Saturation	84.4%
Q30 Bases in Barcode	96.4%
Q30 Bases in RNA Read	59.4%
Q30 Bases in UMI	96.5%

Mapping ?	
Reads Mapped to Genome	95.8%
Reads Mapped Confidently to Genome	92.9%
Reads Mapped Confidently to Intergenic Regions	5.2%
Reads Mapped Confidently to Intronic Regions	25.5%
Reads Mapped Confidently to Exonic Regions	62.2%
Reads Mapped Confidently to Transcriptome	58.2%
Reads Mapped Antisense to Gene	1.2%



Sample	
Sample ID	ETV6-RUNX1_1
Sample Description	
Chemistry	Single Cell 3' v2
Include introns	False
Reference Path	...nger/refdata-cellranger-GRCh38-3.0.0
Transcriptome	GRCh38-3.0.0
Pipeline Version	cellranger-6.0.1

Cell calling




Background GEMs: low #UMI/barcode

Other parameters




Loupe browser



Summary Analysis


t-SNE Projection




Top Features by Cluster (Log2 fold-change, p-value)

Feature	ID	Name	Cluster 1		Cluster 2		Cluster 3		Cluster 4		Cluster 5		Cluster 6		Cluster 7	Cluster 8	Cluster 9
			L2FC	p-value	L2FC	p-value	L2FC	p-value	L2FC	p-value	L2FC	p-value	L2FC	p-value			
ENSG00000277734	TRAC		1.04	1e-1	-0.55	1e+0	-0.20	1e+0	-1.88	4e-1	-0.37	1e+0	-0.72	1e+0	-0.72	1e+0	-0.72
ENSG00000008517	IL32		1.04	9e-1	-5.18	7e-2	-0.54	1e+0	-3.48	1e+0	-3.27	7e-1	-2.31	1e+0	-2.31	1e+0	-2.31
ENSG00000124785	NRN1		1.02	2e-2	0.00	1e+0	-0.08	1e+0	-1.95	3e-2	-1.25	2e-1	-1.49	1e+0	-1.49	1e+0	-1.49
ENSG00000104870	FCGRT		0.92	5e-2	-1.35	6e-2	-0.32	1e+0	-0.40	1e+0	-1.30	2e-1	-1.07	1e+0	-1.07	1e+0	-1.07
ENSG00000185650	ZFP36L1		0.84	1e-1	-0.43	1e+0	-0.51	8e-1	-0.60	1e+0	-0.66	1e+0	-0.48	1e+0	-0.48	1e+0	-0.48
ENSG00000136929	HEMGN		0.82	1e-1	-0.54	1e+0	0.31	9e-1	-0.76	1e+0	-1.30	2e-1	-1.01	1e+0	-1.01	1e+0	-1.01
ENSG00000214113	LYRM4		0.75	2e-1	0.30	1e+0	-0.29	1e+0	-0.89	1e+0	-0.72	1e+0	-1.69	1e+0	-1.69	1e+0	-1.69
ENSG00000101197	BIRC7		0.74	2e-1	0.23	1e+0	-0.81	2e-1	0.45	1e+0	-0.47	1e+0	1.00	8e-1	-1.1	1e+0	-1.1
ENSG00000172889	EGFL7		0.68	4e-1	-0.25	1e+0	-0.19	1e+0	-0.32	1e+0	-0.54	1e+0	-0.44	1e+0	-0.44	1e+0	-0.44
ENSG00000096384	HSP90AB1		0.66	4e-1	-0.35	1e+0	-0.17	1e+0	-0.32	1e+0	-0.03	1e+0	-0.79	1e+0	-0.79	1e+0	-0.79

Sequencing Saturation



Median Genes per Cell



10x single cell portfolio



Flex

Fix, batch, and run on your schedule.

Protein coding gene coverage

Performs well with low-quality and FFPE samples

Multiomic readouts from the same cell

Gene expression
Protein
CRISPR

Throughput options

Up to 8M cells (1–3,072 samples) per run

Products

[Flex](#) >



Universal

Species agnostic. Maximum versatility.

Whole transcriptome coverage

Delivers broadest set of information, including isoforms, SNPs, etc.

Multiomic readouts from the same cell

3' or 5' gene expression
TCR/BCR
Protein
CRISPR

Throughput options

Up to 160,000 cells (1–8 samples) per run

Products

[Universal 3'](#) >

[Universal 5'](#) >



Epi Chromatin

Unmask epigenomic profiles.

ATAC-seq chemistry

Explore open chromatin regions
Link directly to 3' gene expression (Multiome kit)

Multiomic readouts from the same cell

Chromatin accessibility
3' gene expression

Throughput options

Up to 80,000 nuclei (1–8 samples) per run

Products

[Epi ATAC](#) >

[Epi Multiome](#) >

Thanks for your attention!